


Niti-S Biliary Uncovered Stent [M-Type]

for malignant biliary strictures



Features

- **Braided and a weaving construction**
- **Radiopaque marker:** 3 (three) at both ends & 2 (two) in the middle

Ordering Information

Percutaneous							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
TN0604	6	4	7	50			
TN0605		5					
TN0606		6					
TN0607		7					
TN0608		8					
TN0609		9					
TN0610		10					
TN0612		12					
TN0804		8			4	7	50
TN0805					5		
TN0806	6						
TN0807	7						
TN0808	8						
TN0809	9						
TN0810	10						
TN0812	12						
TN1004	10		4	7	50		
TN1005			5				
TN1006		6					
TN1007		7					
TN1008		8					
TN1009		9					
TN1010		10					
TN1012		12					

Percutaneous							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
TN0604-6	6	4	6	50			
TN0605-6		5					
TN0606-6		6					
TN0607-6		7					
TN0608-6		8					
TN0609-6		9					
TN0610-6		10					
TN0612-6		12					
TN0804-6		8			4	6	50
TN0805-6					5		
TN0806-6	6						
TN0807-6	7						
TN0808-6	8						
TN0809-6	9						
TN0810-6	10						
TN0812-6	12						

Ordering Information

Endoscopic							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
BN0604	6	4	7	180			
BN0605		5					
BN0606		6					
BN0607		7					
BN0608		8					
BN0609		9					
BN0610		10					
BN0612		12					
BN0804		8			4	7	180
BN0805					5		
BN0806	6						
BN0807	7						
BN0808	8						
BN0809	9						
BN0810	10						
BN0812	12						
BN1004	10		4	7	180		
BN1005			5				
BN1006		6					
BN1007		7					
BN1008		8					
BN1009		9					
BN1010		10					
BN1012		12					

Endoscopic							
Code	Stent		Delivery				
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)			
BN0604-6	6	4	6	180			
BN0605-6		5					
BN0606-6		6					
BN0607-6		7					
BN0608-6		8					
BN0609-6		9					
BN0610-6		10					
BN0612-6		12					
BN0804-6		8			4	6	180
BN0805-6					5		
BN0806-6	6						
BN0807-6	7						
BN0808-6	8						
BN0809-6	9						
BN0810-6	10						
BN0812-6	12						